


Product Catalogue for Life Science

Cell counters & Live cell imaging systems





NanoEntek

www.nanoentek.com
sales@nanoentek.com

NanoEntek, Inc. (Head Office)

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel: +82-2-6220-7940 / Fax: +82-2-6220-7999

NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA
Tel:+1-781-472-2558 / Fax: +1-781-790-5649

Contents

About NanoEntek

Company Overview	4~5
Partnership & Global Networks	6

Product Lineup	8
-----------------------------	---

Cell Counter

ADAM™ MC2	10
ADAM™ MC Plus	11
ADAM™ CellIT	12
ADAM™ CellIT Plus	13
ADAMII™ LS	14
ADAMII™ CDx	15
ADAM™ SCC2	16
EVE™	17
EVE™ PLUS	18
EVE™ HT	19
EVE™ HT FL	20

Live Cell Imaging System

JuLI™ Stage	21
-------------------	----

Disposable Hemocytometer

C-Chip 2ch	22
C-Chip 4ch	23
S-Chip	24

Ordering Information	26~28
-----------------------------------	-------

About NanoEntek

Outline

Company Overview
Partnership & Global Networks

Company Overview




Digital Bio-It all began with young scientists willing to contribute to society

In 2000, NanoEntek(KOSDAQ: 039860) was established initially as Digital Bio with 10 scientists who were in the field of various sciences in Seoul National University. NanoEntek possesses the technology for a key platform called “Lab-On-a-Chip”, which integrates laboratory processes on a microfluidic chip with the size of a fingernail, making significant contributions to the fields of medical and life science engineering.

Our mission : Support researcher for a better world

From the beginning until today, NanoEntek is aiming to develop innovative and easy to use products (automated cell counter and analyzer) for life science and cell therapy.

We help them consume their valuable their time and passion for discovery, instead of repetitive and tedious work. NanoEntek always keeps our initial mission in mind and help researchers contribute for the better world.





Success of cell counters

The first successful products was Countess™, the world's first personal cell counter which was a shift in paradigm of cell counters. Following this success, it developed the series of other cell counters; EVE series for academic laboratory, ADAM™ CellT cell counter for cGMP facility, and ADAMII™ LS novel fluorescence cell counter/analyzer.

The second success came from ADAM™ rWBC, the world's first benchtop residual white blood cell counter with FDA 510(k) cleared. It has been selected as the QC standard equipment of the American Red Cross (ARC) and installed in more than 200 blood banks globally.

Through the strong performance (high-accuracy, speed and usability) and customer's satisfaction, we launched 6 new products in 2023 including EVE™ HT FL (High-Throughput Fluorescence Cell Counter) and JuLI™ Stage (Real-time live cell imaging system). They can even fulfill all your needs as well from the research level to the large-scale multi-site manufacturing.

2009

- EVE™ is launched
- Patents is sold to Life Technologies (USD 13M)
- Contract new products development with Life Technologies

2010

- Receives the "USD 5 Million Export Prize" presented by Korea International Trade Association
- Receives the "IR52 JangYoungShil Award" for Countess

2011

- Strategic partnership with SK telecom (invest US 13M\$)

2013

- JuLI™ Br is launched

2016

- SK telecom become the largest share holder
- JuLI™ Stage is launched

Partnership & Global Networks

Accelerate an effort to secure a global sales network to hold partnership with world's largest biomedical device production and distribution company 'ThermoFisher Scientific' and to secure a global distribution network with 150 nations.



- 2008 : Countess supply contract
- 2009 : Microporator patent sold
- 2009 : New product development supply contract



- 2019 : Partnership & supply contract
- 2020 : Supply Millicell® Disposable Hemocytometer

Global Market

Secure sales network through partnership with international global compnies

Hold distribution network with 150 nations

Global Partners

Our worldwide network with global partners



Our Valuable Customers

Our valuable customers who trust and support us for a long time



Product Lineup

Life Science Research

Cell Counter

Live Cell Imaging System

Disposable Hemocytometer

Life Science Research

Cell Counter

NEW



ADAM™ MC Plus

NEW



ADAM™ CellIT Plus



ADAM™ MC2



ADAM™ CellIT

NEW



EVE™ HT FL

NEW



EVE™ HT



EVE™



EVE™ PLUS

NEW



ADAMII™ LS

NEW



ADAMII™ CDx



ADAM™ SCC2

Live Cell Imaging System

NEW



JuLI™ Stage

Disposable Hemocytometer



C-Chip 2ch



C-Chip 4ch



S-Chip

Life Science Research

Outline

Cell Counter

Live Cell Imaging System

Disposable Hemocytometer

ADAM™ MC2 PI/AO staining method



Accurate Automated Cell Counter

ADAM™ MC2 is an accurate automated fluorescence cell counter that performs viability and cell counting measurements within 25 seconds.

KEY FEATURES & BENEFITS

- Less than 25 seconds to get results
- Higher accuracy
- Sensitive CMOS detection
- WiFi available

Applications

- Cell counting
- Viability test
- Cell Therapeutic Applications



Specifications



AccuChip Kit (PI & AO/PI)

Analysis time	< 25 sec / test	Measuring range	5 X 10 ⁴ ~ 4 X 10 ⁷ cells/mL (PI) 5 X 10 ⁴ ~ 2 X 10 ⁷ cells/mL (AO/PI)
Focus	Auto	Loading sample volume	13 µL / test (AccuChip 4X)
Cell type	Cell line(clumpy cell, single cell), PBMCs, Adipose Stem cell, Total white blood cell, Primary cells	Measuring sample volume	3.4 µL / test (AccuChip 4X)
Weight	7.0 kg	Staining solution	Propidium Iodide stain Acridine Orange stain
Dimensions	277 x 276 x 270 mm (W x D x H)	Connectivity	WiFi (USB Dongle)



ADAM™ MC Plus ●● 2 FLs (AO/DAPI) & Bright field

Most Accurate Fluorescence Cell Counter

ADAM™ MC Plus is the new standard of automated fluorescence cell counter. It measures the number of total cells, viable cells, non-viable cells and shows viability results. In addition, it analyzes the cell size and aggregation ratio as well.

KEY FEATURES & BENEFITS

- Only 15 μ L of sample volume needed
- 4 tests available with 1 slide
- 45 sec per test
- Low CV value
- PBMCs, Stem Cells, Primary Cells, Cell Lines

Applications

- Cell counting
- Viability test
- Cell therapeutic applications
- Cell size
- Cell aggregation ratio



AccuPlus Slide & Reagent


Specifications

Analysis time	10 sec / test (fast mode) 30 sec / test (cell size mode)
Focus	Auto
Cell type	Cell line(clumpy cell, single cell), PBMCs, Adipose Stem cell, Total white blood cell, Primary cells
Weight	7.0 kg
Dimensions	277 x 276 x 270 mm (W x D x H)

Measuring range	5 X 10 ⁴ ~ 2 X 10 ⁷ cells/mL (Detectable) 4 X 10 ⁵ ~ 1 X 10 ⁷ cells/mL (Optimal)
Loading sample volume	15 μ L / test
Measuring sample volume	3.2 μ L / test
Staining solution	Acridine orange / DAPI solution

ADAM™ CellIT

PI/AO staining method



Automated Cell Counter for cGMP

ADAM™ CellIT is an automated cell counter that is available in cGMP production environment.

ADAM™ CellIT complies with 21 CFR Part 11 which is a regulation on electronic records and signatures for use in cGMP facilities. The data cannot be modified by any user.

KEY FEATURES & BENEFITS

- 21 CFR Part 11 Compliance
- Electronic signatures are required for use in cGMP facilities
- WiFi available

Applications

- Cell counting
- Viability test
- Cell Therapeutic Applications



Specifications



AccuChip Kit (PI & AO/PI)

Analysis time	< 25 sec / test	Measuring range	5 X 10 ⁴ ~ 4 X 10 ⁶ cells/mL (PI) 5 X 10 ⁴ ~ 2 X 10 ⁷ cells/mL (AO/PI)
Focus	Auto	Loading sample volume	13 µL / test (AccuChip 4X)
Cell type	Cell line(clumpy cell, single cell), PBMCs, Adipose Stem cell, Total white blood cell, Primary cells	Measuring sample volume	3.4 µL / test (AccuChip 4X)
Weight	7.0 kg	Staining solution	Propidium Iodide stain Acridine Orange stain
Dimensions	277 x 276 x 270 mm (W x D x H)	Connectivity	WiFi (USB Dongle)



ADAM™ CellIT Plus ●● 2 FLs (AO/DAPI) & Bright field

Most Accurate Fluorescence Cell Counter for cGMP

ADAM™ CellIT Plus is the new standard of automated fluorescence cell counter. It measures the number of total cells, viable cells, non-viable cells and shows viability results. In addition, it analyzes the cell size and aggregation ratio as well.

KEY FEATURES & BENEFITS

- Only 15 μ L of sample volume needed
- 4 tests available with 1 slide
- 45 sec per test
- Low CV value
- PBMCs, Adipocytes, Leukocytes, Stem Cells, Hepatocytes, Splenocytes, Cell Lines

Applications

- Cell counting
- Viability test
- Cell therapeutic applications
- Cell size
- Cell aggregation ratio

NEW



AccuPlus Slide & Reagent

Specifications

Analysis time	10 sec / test (Fast mode) 30 sec / test (Cell size mode)
Focus	Auto
Cell type	Cell line(clumpy cell, single cell), PBMCs, Adipose Stem cell, Total white blood cell, Primary cells
Weight	7.0 kg
Dimensions	277 x 276 x 270 mm (W x D x H)

Measuring range	$5 \times 10^4 \sim 2 \times 10^7$ cells/mL (Detectable) $4 \times 10^5 \sim 1 \times 10^7$ cells/mL (Optimal)
Loading sample volume	15 μ L / test
Measuring sample volume	3.2 μ L / test
Staining solution	Acridine Orange stain DAPI stain

ADAMII™ LS

● ● ● 3 FLs (GFP/RFP/DAPI) & Bright field



Image-based fluorescence cell analyzer

ADAMII™ LS is a benchtop fluorescent intensity based cell analyzer providing versatile assays and high accurate results (similarity to flow cytometer).

KEY FEATURES & BENEFITS

- Fluorescence level
- Transfection efficiency
- Protein expression
- Dot plot
- 21 CFR part 11 compliance for GMP

Applications

- Fluorescence expression
- Cell cycle
- Apoptosis
- Total cell counting & Viability
- Cell size

NEW



Assay Slide & Reagents

Specifications

Analysis time	App. 2 min ~ 4 min 30 sec
Objective lens	10 X
Operating power	100 – 240 V, 1.5 A, 50/60 Hz
Weight	19.3 kg
Dimensions	300 x 420 x 370 mm (W x D x H)

Measurement range	5 X 10 ⁴ ~ 4 X 10 ⁶ cells/mL (PI) 5 X 10 ⁴ ~ 2 X 10 ⁷ cells/mL (AO/PI)
Loading volume	25 µL
Measuring volume	≤ 7.8 µL
Electronic input	12VDC, 5.0A
Light source	Bright field, UV, Blue, Green LED



ADAMII™ CDx ● ● ● 3 FLs (UV/Blue/Green) & Bright field

Image-based immune cell phenotyping instrument

ADAMII™ CDx, all-in-one system with 4-channel (Bright field and 3 fluorescent channels) for cell therapy R&D and manufacturing process.

KEY FEATURES & BENEFITS

- Absolute CD + cell counting
- Various CD marker panels
- Total cell counting, viability, cell size, growth curve
- Apoptosis assay
- 21 CFR part 11 compliance for GMP

Applications

Immune cell phenotyping

- T cell : CD34/CD8/CD3
- NK cell : CD3/CD56&16/CD45
- B cell : CD3/CD19/CD45
- Monocyte : CD3/CD14/CD45

NEW



Specifications

Analysis time	App. 2 min ~ 4 min 30 sec	Measurement range	$5 \times 10^4 \sim 5 \times 10^6$ cells/mL
Objective lens	10 X	Loading volume	25 μ L
Weight	19.3 kg	Measuring volume	$\leq 7.8 \mu$ L
Dimensions	300 x 420 x 370 mm (W x D x H)	Light source	Bright field, UV, Blue, Green LED

ADAM™ SCC2

PI staining method



Fast Automated Somatic Cell Counter

ADAM™ SCC2 is a fast, easy and accurate automated somatic cell counter that counts somatic cell to determine raw milk quality.

KEY FEATURES & BENEFITS

- Less than 13 seconds for a test (4ch)
- Optimized image analysis software
- WiFi available

Applications

- Counting somatic cells
- Raw milk quality control



Specifications

Analysis time	< 13 sec / test	Measuring range	0.05 ~ 1.15 x 10 ⁶ cells/mL
Focus	Auto	Loading sample volume	13 µL / test (SomaChip 4X)
Cell type	Somatic cell	Measuring sample volume	3.4 µL / test (SomaChip 4X)
Weight	7.0 kg	Staining solution	Propidium Iodide stain
Dimensions	277 x 276 x 270 mm (W x D x H)	Connectivity	WiFi (USB Dongle)



Somatic Cell Count Kit

EVE™  TB staining method



Automated Cell Counter

EVE™ uses state-of-the-art optics and image analysis to automate cell counting.

The EVE™ is a benchtop counter, designed to measure cell count and viability (live, dead, and total cells) accurately and precisely, using the standard trypan blue technique.

KEY FEATURES & BENEFITS

- Fine distinction of clumped cell
- Cell size gating
- Data storage & analysis
- No need of maintenance
- Bench top size & LCD touch screen
- Manual focus

Applications

- Cell viability (%)
- Total cell counting
- Cell concentration measurement
- Data analysis with graphical reports



EVE™ Slide

Specifications

Analysis time	< 20 sec / test
Power	100-240V~, 50/60Hz
Display	7" LCD touch screen
Data storage	USB drive (4 GB)
Data format	JPEG (Image) , CSV (Raw data)
Weight	2.1 kg
Dimensions	270 x 200 x 190 mm (W x D x H)

Measuring range	1 X 10 ⁴ ~ 1 X 10 ⁷ cells/mL (Detectable) 1 X 10 ⁵ ~ 4 X 10 ⁶ cells/mL (Optimal)
Loading sample volume	10 µL / test
Cell size range	5 ~ 60 µm
Staining solution	Trypan blue stain

EVE™ PLUS TB staining method



The World's Fastest Automated Cell Counter

EVE™ Plus identifies and counts clumpy cells individually for accurate analysis within 1 second.

KEY FEATURES & BENEFITS

- Less than 1 second to get results
- Easy to use
- Automatically saved up to 500 results
- Fine distinction of clumped cells
- 2 focus mode: auto, manual
- WiFi available



Applications

- Total cell counting (live, dead)
- Viability (%)
- Cell size

Specifications



EVE™ Slide

Analysis time	<1 sec (manual) / <10 sec (auto)	Measuring range	1x10 ⁴ ~ 2x10 ⁷ cells/mL
Focus	Auto & Manual	Loading sample volume	10 µL / test
Connectivity	WiFi (USB Dongle)	Cell size range	5 ~ 60 µm
Cell type	Cell line (clumpy cell, single cell)	Staining solution	Trypan blue stain
Weight	4.0 kg		
Dimensions	274 x 274 x 333 mm (W x D x H)		

EVE™ HT

TB staining method



High-throughput Automated Cell Counter

The EVE™ HT is a high-throughput automated multiple cell counter, providing 48-sample counting in just 2.5 minutes to analyze primary cells and cell lines. With small sample volume and extensive capacity, the EVE™ HT provides accurate cell counting and viability measurements.

KEY FEATURES & BENEFITS

- Up to 48 samples per count
- 3 minutes to count 48 samples
- 20 µL required for each sample
- User-independent consistency in results
- Cell declustering technique
- 21 CFR Part 11 Compliance

Applications

- Total cell counting (viable, non-viable)
- Viability (%)
- Cell size

NEW



Specifications



EVE™ HT Counting Kit

Analysis time	< 3 min / 48 tests	Measuring range	$1 \times 10^4 \sim 1 \times 10^7$ cells/mL
Power	100-240V~, 50/60Hz	Cell size range	5 ~ 80 µm
Weight	58 kg	Loading sample volume	20 µL / channel
Dimensions	586 x 477 x 458 mm (W x L x H)	Staining solution	Trypan blue stain

EVE™ HT FL

●● 2 FLs (AO/DAPI) & Bright field




High-throughput Fluorescence Cell Counter

The EVE™ HT FL is a high-throughput automated fluorescence cell counter equipped with bright field and two fluorescence channels (AO/DAPI). In just 3 minutes, up to 48 samples can be counted and analyzed. EVE™ HT FL delivers brilliant precision and accuracy, making it the optimal choice for both cell lines and primary cell analysis in a variety of applications.

KEY FEATURES & BENEFITS

- Dual-Fluorescence (AO/DAPI) & Bright Field
- 48 Samples in just 3 minutes
- Only 20 µL sample volume required
- Real Cell Size & Viability
- Accurately identify and count clumpy cells

NEW



EVE™ HT FL Counting Kit

Specifications

Analysis time	3 ~ 6 minutes for 48 samples	Measuring range	Detectable : $1 \times 10^4 \sim 2 \times 10^7$ cells/mL Optimal : $1 \times 10^5 \sim 1 \times 10^7$ cells/mL
Power	100 ~ 240V, 50/60Hz	Loading sample volume	20 µL / channel
Cell size range	5 ~ 85 µm (Detectable) 5 ~ 80 µm (Optimal)	Staining solution	AO/DAPI mixed solution
Weight	61 kg	Channel	Dual fluorescence channels (AO&DAPI)
Dimensions	586 x 461 x 458 mm (W x D x H)	Operation system	Windows 10
		21 CFR Part 11 compliance	Available (Optional)

JuLI™ Stage 3 FLs (GFP/RFP/DAPI)




Real-time Live Cell Imaging System

JuLI™ Stage, the new standard of real-time live cell imaging system is designed to get time-lapse images and make taking cell videos much easier.

KEY FEATURES & BENEFITS

- Compact and compatible with a standard CO₂ incubator
- Fully automated X-Y-Z stage
- Multi-channel fluorescence imaging (GFP, RFP, DAPI and Bright field)
- Easy & powerful software
- Take and analyze image in real-time

Applications



Specifications

Light source	Blue, Green, UV LED (Intensity adjustable)
Fluorescence	DAPI, GFP, RFP
Stage	Automated, motorized, X-Y-Z stage / Interchangeable vessel holder
PC	TBD
Electronic input	12VDC, 5.0A
Weight	18.5 kg / 41 lbs
Dimensions	429 x 310 x 324 mm (W x D x H)

Objective lens	4X, 10X, 20X + Digital zoom Inter-changeable objective lens
Camera	High-sensitivity monochrome CCD (Sony sensor 2/3") / 1,936 x 1,456 pixels (2.8M), 53 FPS, 14 bit
Exported formats	Image : JPEG, TIFF, BMP, PNG Video : AVI Raw data : CSV
Operating power	100-240V, 1.5A, 50/60Hz
Operating environment	5-40°C, 20-95%

C-Chip 2ch



The World's First True Disposable Hemocytometer

A precision disposable plastic hemocytometer, developed to solve the problems of conventional glass hemocytometer.

KEY FEATURES & BENEFITS

- 2 tests in 1 slide
- No coverslip required
- No more exposure to hazardous materials (e.g., AIDS patient blood, urine)
- Accurate and reliable
- Quartz grade optical plastic



Applications

- Blood analysis
- Cell culture
- Microbiology

Specifications



Cat. No.	DHC-N01	DHC-B01	DHC-B02	DHC-F01	DHC-M01
Grid	Neubauer Improved	Bürker	Bürker-Türk	Fuchs Rosenthal	Malassez
Unit	50 Slide	50 Slide	50 Slide	50 Slide	50 Slide
Chamber depth	100 µm	100 µm	100 µm	200 µm	200 µm
Loading volume	10 µL	10 µL	10 µL	20 µL	20 µL
Dimensions	25 x 75 x 1.6 mm (W x D x H)				

C-Chip 4ch



Disposable Hemocytometer 4 channel

4 channel precision disposable plastic hemocytometer, developed to solve the problem of conventional glass hemocytometer.

KEY FEATURES & BENEFITS

- 4 tests in 1 slide
- No coverslip required
- No more exposure to hazardous materials
- Accurate and reliable
- Quartz grade optical plastic

Applications

- Blood analysis
- Cell culture
- Microbiology



Specifications

Cat. No.	DHC-N04	DHC-B04	DHC-F04
Grid	Neubauer Improved	Bürker	Fuchs Rosenthal
Unit	50 Slide	50 Slide	50 Slide
Chamber depth	100 μ m	100 μ m	200 μ m
Loading volume	10 μ L	10 μ L	20 μ L
Dimensions	25 x 75 x 1.6 mm (W x D x H)		

S-Chip

CE

Disposable Counting Chamber for Sperm

A precise and disposable plastic chamber with two ports for sperm counting.

KEY FEATURES & BENEFITS

- 2 tests in 1 slide
- No coverslip required
- No need to wash for reuse
- Accurate and reliable
- Quartz grade optical plastic

Applications

- IVF, IUI: Sperm counting



Specifications

Chamber depth	20 μ m
Loading volume	2 ~ 3 μ L
Capacity	50 slides / box (100 counts)
Dimensions	25 x 75 x 2 mm (W x D x H)

Ordering Information

Outline

Cell Counter

Live Cell Imaging System

Disposable Hemocytometer

Cell Counter

ADAM™ Series

Catalog Number	Product Name	Specification
ADAM-MC2	Automated cell counter, ADAM™ MC2	Main device
ADAM-CellT	Automated cell counter for cGMP, ADAM™ CellT	Main device
AD4K-200	AccuChip kit (PI staining)	- AccuChip 4X x 200 pcs - T(PI) sol. x 2 ea, N sol. x 1 ea
AD4K-200AO	AccuChip kit (AO/PI staining)	- AccuChip 4X x 200 pcs - T(AO) sol. x 2 ea, N sol. x 1 ea
ADR-1000	Accustain solution (PI staining)	T(PI) sol. x 4 ea, N sol. x 2 ea
ADR-1000AO	Accustain solution (AO/PI staining)	T(AO) sol. x 4 ea, N sol. x 2 ea
ADB-500	ADAM calibration beads	Calibration beads 5 mL x 1 bottle
QCS-001	QC Slide for ADAM™ CellT	
ADAM-MC Plus	Automated fluorescence cell counter, ADAM™ MC Plus	Main device
ADAM-CellT Plus	Automated fluorescence cell counter for cGMP, ADAM™ CellT Plus	- Main device - 21 CFR PART 11 appendix
APAD-400	Cell viability reagent	1 mL x 6 tubes (400 tests)
APAS-100	AccuPlus Slide 4ch	4ch. Slide 100 ea
ADAM-SCC2	Automated somatic cell counter, ADAM™ SCC2	
CRS-K02	Somatic cell count kit	- 4ch chip x 100 slides (400 tests) - SCC sol. x 2 ea

ADAMII™ Series

Catalog Number	Product Name	Specification
ADAMII-LS	Image-based fluorescence cell analyzer, ADAMII™ LS	Main device, Laptop
ALAD-100	Cell viability reagent (AO/DAPI staining)	AO/DAPI stain 0.5 mL x 2 tubes (100 tests)
ALPI-100	PI cell cycle reagent	PI stain 1.25 mL x 2 tubes (100 tests)
ALAP-100	Apoptosis detection kit	- AnnexinV-PE stain: 0.5 mL x 1 tube (100 Tests) - DAPI solution: 125 µL x 1 tube (100 Tests) - AnnexinV binding buffer: 10 mL x 1 tube (100 tests)
A2AS-051	ADAMII Assay slide	1 ch x 50 slides/ case
ADAMII-CDx	Image-based immunophenotyping instrument, ADAMII™ CDx	Main device, Laptop
ATK-025	ADAMII™ T cell assay kit(CD4/CD8/CD3)	- 25 tests / T reagent, 10X RBC lysis buffer, Calibration beads, Assay slide
ABK-025	ADAMII™ B cell assay kit (CD3/CD19/CD45)	- 25 tests / B reagent, 10X RBC lysis buffer, Calibration beads, Assay slide
ANK-025	ADAMII™ NK cell assay kit (CD3/CD16+CD56/CD45)	- 25 tests / NK reagent, 10X RBC lysis buffer, Calibration beads, Assay slide
AMK-025	ADAMII™ Monocyte assay kit (CD3/CD14/CD45)	- 25 tests / Monocyte reagent, 10X RBC lysis buffer, Calibration beads, Assay slide
ACVK-025	ADAMII™ Cell viability kit	- 25 tests / Cell viability reagent, Assay slide

EVE™ Series

Catalog Number	Product Name	Specification
EVE-MC	Automated cell counter, EVE™	Main device, Desktop computer
EVE-MC2	The World's fastest cell counter, EVE™ Plus	Main device, Desktop computer
EVS-050	EVE cell counting slide	50 slides (100 tests)
EBB-001	Test beads	1 mL vial x 1 ea
EVE HT	High-throughput automated cell counter, EVE™ HT	Main device, Desktop computer
EVH-020	EVET™ HT counting kit	- 960 tests / kit - Counting plate, Mixing well plate, Trypan blue stain 0.4 %, Reservoir
EHPQ-001	EVE™ HT QC plate, Low level	
EHPQ-002	EVE™ HT QC plate, Middle level	
EHPQ-003	EVE™ HT QC plate, High level	
EVE HT FL	High-throughput fluorescence cell counter, EVE™ HT FL	Main device, Desktop computer
EVFL-020	EVE™ HT FL Counting kit	- 960 tests / kit - Counting plate, Mixing well plate, Reservoir
EVAD-960	AO/DAPI Staining solution	AO/DAPI 20 mL x 2 bottles

Live Cell Imaging System

JuLI™ Series

Catalog Number	Product Name	Specification
JS1000S	JuLI™ Stage, Starter pack	JuLI Stage basic set (JS1000), Desktop computer (JP0200), 3 Objective lenses (4X, 10X & 20X)
JS1000	JuLI™ Stage, Real-time live cell imaging system	Main Device, power supply, control box
JP0100	Desktop computer	- CPU: Intel i5, 9 generation or over spec. - OS: Windows 10 Pro 64 bit - RAM: 16 GB - Hard drive: 2 TB - Network: Gigabit Ethernet
JMO100	Desktop computer	24" Full HD (1920 x 1980) monitor
JP0150	External hard disk drive (Optional)	Total 8 TB (4 TB x 2 ea)
JO0004	Objective lens (4X)	Magnification : 4X, NA : 0.16
JO0010	Objective lens (10X)	Magnification : 10X, NA : 0.3
JO0020	Objective lens (20X)	Magnification : 20 X, NA : 0.45
JVH001	Vessel holder (Optional)	Micro Slide (26 x 76 mm)
JVH002	Vessel holder (Optional)	Petri Dish (35 mm)
JVH003	Vessel holder (Optional)	Petri Dish (60 mm)
JVH004	Vessel holder (Optional)	Petri Dish (100 mm)
JVH005	Vessel holder (Optional)	T-Flask (25 & 75cm ²)
JSCT100	JuLI analysis software (Scratch)	- JuLI Scratch STAT - JuLI Scratcher
JSPT100	JuLI analysis software (Spheroid)	JuLI Spheroid STAT

Disposable Hemocytometer

C-Chip 2ch

Catalog Number	Product Name	Contents
DHC-N01	Neubauer Improved, C-chip	100 μm (depth), 50 slides
DHC-F01	Fuchs-Rosenthal, C-chip	200 μm (depth), 50 slides
DHC-B01	Bürker, C-chip	100 μm (depth), 50 slides
DHC-B02	Bürker-Türk, C-chip	100 μm (depth), 50 slides
DHC-M01	Malassez, C-chip	200 μm (depth), 50 slides

C-Chip 4ch

Catalog Number	Product Name	Contents
DHC-N04	Neubauer Improved, C-chip 4ch	100 μm (depth), 50 slides
DHC-B04	Bürker-Türk, C-chip 4ch	100 μm (depth), 50 slides
DHC-F04	Fuchs-Rosenthal, C-chip 4ch	200 μm (depth), 50 slides

S-Chip

Catalog Number	Product Name	Contents
DCS-S01	Disposable plastic chamber for sperm counting, S-chip	20 μm (depth), 50 slides

MEMO



www.nanoentek.com
sales@nanoentek.com

NanoEntek, Inc. (Head Office)

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel: +82-2-6220-7940 / Fax: +82-2-6220-7999

NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA
Tel:+1-781-472-2558 / Fax: +1-781-790-5649



SCAN ME!

Homepage www.nanoentek.com

Facebook facebook.com/nanoentek

Instagram instagram.com/nanoentek

LinkedIn linkedin.com/company/nanoentek

Youtube www.youtube.com/c/NanoEntek

Blog www.blog-nanoentek.com



e-mail sales@nanoentek.com
website www.nanoentek.com

NanoEntek, Inc. (Head Office)

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel: +82-2-6220-7940 / Fax: +82-2-6220-7999

NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA
Tel: +1-781-472-2558 / Fax: +1-781-790-5649